



SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE ADULT FINGER CLIP FOR MP3



Réf. 222260



SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE ADULT FINGER CLIP FOR MP3



Réf. **222260**

FEATURES			
Medical device	No		
Marquage	Barcode	3700446030131	
Guarantee (month):	12	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,105
SPECIFIC STORAGE CONDITION			
Temperature min	1	Temperature max	50
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:41 - This article is available on our website : www.gsh-med.fr with the following reference 222260

All characteristics are given for information only and may be changed