

SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE PEDIATRIC FINGER CLIP FOR MP3

Réf. 222261



SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE PEDIATRIC FINGER CLIP FOR MP3



Réf. 222261

FEATURES			
Medical device	No		
Marquage		Barcode	3700446030148
Guarantee (month):	12	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,094
SPECIFIC STORAGE CONDITION			
Temperature min	1	Temperature max	50
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:41 - This article is available on our website : www.gsh-med.fr with the following reference 222261

All characteristics are given for information only and may be changed