

CARRYING BAG FOR OXIMETER MASTER PALM 2&3

Réf. 223003



CE
0123



CARRYING BAG FOR OXIMETER MASTER PALM 2&3



Réf. **223003**

FEATURES				
Medical device	Yes	Classe IIb	Accuracy measurement	Classe II
Marquage	CE0123		Barcode	3700446025441
Guarantee (month):	0		Customs Code	90330000
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :	0,035	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



31/03/2022 01:41 - This article is available on our website : www.gsh-med.fr with the following reference 223003

All characteristics are given for information only and may be changed