

SILICON PROTECT COVER FOR MASTER PALM 2/3 (2 PCS)

Réf. 223004



CE
0123



SILICON PROTECT COVER FOR MASTER PALM 2/3 (2 PCS)



Réf. **223004**

FEATURES				
Medical device	Yes	Classe IIb	Accuracy measurement	Classe II
Marquage	CE0123		Barcode	3700446025434
Guarantee (month):	0		Customs Code	90330000
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :	0,074	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)		Temperature max	(N/A)
Humidity min	(N/A)		Humidity max	(N/A)



31/03/2022 01:41 - This article is available on our website : www.gsh-med.fr with the following reference 223004

All characteristics are given for information only and may be changed