


## CARDIOMATE 12 - ELECTROCARDIOGRAPH 12 CHANNELS

Réf. 260012



# CARDIOMATE 12 - ELECTROCARDIOGRAPH 12 CHANNELS



**Réf. 260012**

ECG Cardiomate 12® are robust and reliable ECG. Featuring a 10" color screen tilt they allow viewing 12 lead simultaneously. Simple to use, recording and printing ECG are triggered with only one button. The Cardiomate 12® are Cardiomate equipped with the latest technology storage and communication (SD, USB card, network) that can record internally to 10 000 ECG, send them to a computer (optional software) or directly to an external printer. The Cardiomate 12® are equipped with automatic interpretation and pacemaker detection. The internal battery life of 3 hours, their lightness and carrying handle make Cardiomate 12® the perfect companion in the office or traveling.- Carrying handheld.- Tilting screen :10-inch color LCD display (resolution: 800x480 px).- SD card included :4GB SD card, 10,000 reports of ECG memory in XML / PDF / BMP / JPG / DICOM.- Autonomy :The operating life is 1080 ECG.- Automatic interpretation and intervals calculation.- Pacemaker detection.- Network connection, USB & printer.- Printing and rapid analysis.- Acquisitions and display 12 simultaneous tracks.- Alphabetic keypad.

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446025496
Guarantee (month):	24	Customs Code	9018110000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	40,0	Width packaging (cm):	
Length (cm) :	42,0	Length packaging (cm) :	
Height (cm) :	9,5	Height packaging (cm) :	
Weight (kg) :	5,200	Weight packaging (kg) :	8,300
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	50
Humidity min	25	Humidity max	85



18/05/2022 21:07 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 260012

All characteristics are given for information only and may be changed