

BAG OF 10 UNITS ECG ADAPTER PRESSURE BUTTON

Réf. 260023



BAG OF 10 UNITS ECG ADAPTER PRESSURE BUTTON



Réf. 260023

FEATURES			
Medical device	No		
Marquage		Barcode	3700446030704
Guarantee (month):	24	Customs Code	90330000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,060
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:42 - This article is available on our website : www.gsh-med.fr with the following reference 260023

All characteristics are given for information only and may be changed