

## 12-CHANNEL ECG PAPER, ROLLED



Réf. 260026



## 12-CHANNEL ECG PAPER, ROLLED



Réf. 260026

FEATURES			
Medical device	No		
Marquage	Barcode	3700446026691	
Guarantee (month):	0	Customs Code	4823400000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	6,0
Length (cm) :		Length packaging (cm) :	22,0
Height (cm) :		Height packaging (cm) :	6,0
Weight (kg) :		Weight packaging (kg) :	0,372
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:42 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 260026

All characteristics are given for information only and may be changed