

MAXI +3 WALL MOUNTED WITH SET OF 3 CUFFS (S, M, L)

Réf. 522707



 Spengler

GARANTIE
2
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE

CLASSE
I_m

MAXI +3 WALL MOUNTED WITH SET OF 3 CUFFS (S, M, L)



Réf. 522707

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0123	Barcode	3700446031152
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	1,154
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/05/2022 12:31 - This article is available on our website : www.gsh-med.fr with the following reference 522707

All characteristics are given for information only and may be changed