

## WALL MOUNTED WITH ADULT CUFF (M)

Réf. 522700



## WALL MOUNTED WITH ADULT CUFF (M)



**Réf. 522700**

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0123	Barcode	3700446028145
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,926
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/05/2022 12:31 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 522700

All characteristics are given for information only and may be changed