

AC ADAPTER FOR AUTOTENSIO TENSONIC BPM

Réf. 527184



AC ADAPTER FOR AUTOTENSIO TENSIONIC BPM



Réf. **527184**

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage	CE0197		Barcode	3700446022839
Guarantee (month):	60		Customs Code	90330000
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :	0,081	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



31/03/2022 01:20 - This article is available on our website : www.gsh-med.fr with the following reference 527184

All characteristics are given for information only and may be changed