

## GREY SPENGLER (S) CUFF (16-24CM) FOR TENSONIC & AUTOTENSIO

Réf. 527462



# GREY SPENGLER (S) CUFF (16-24CM) FOR TENSIONIC & AUTOTENSIO



Réf. **527462**

Accessory for ES-60® AUTOTENSIO® TENSIONIC® & Effitens®

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage	CE0197		Barcode	3700446022969
Guarantee (month):	6		Customs Code	90262040
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :	0,102	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



26/04/2022 12:32 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 527462

All characteristics are given for information only and may be changed