

## STATION DE DIAGNOSTIC SPOT-CHECK VS PNI

Réf. 546000



## STATION DE DIAGNOSTIC SPOT-CHECK VS PNI



Réf. **546000**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446035822
Guarantee (month):	12	Customs Code	9018199560
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	1,489
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:51 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 546000

All characteristics are given for information only and may be changed