

SPO2 SENSOR PEDIATRIC FOR SPOT-CHECK 2VS



Réf. 546011



SPO2 SENSOR PEDIATRIC FOR SPOT-CHECK 2VS



Réf. 546011

FEATURES			
Medical device	No		
Marquage		Barcode	3700446029616
Guarantee (month):	6	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,056
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:51 - This article is available on our website : www.gsh-med.fr with the following reference 546011

All characteristics are given for information only and may be changed