

## STATION DE DIAGNOSTIC SPOT-CHECK 2VS PNI AVEC SOCLE A ROULETTES

Réf. 546002



# STATION DE DIAGNOSTIC SPOT-CHECK 2VS PNI AVEC SOCLE A ROULETTES



Réf. **546002**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446035839
Guarantee (month):	12	Customs Code	9018199560
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	60
Humidity min	10	Humidity max	95



31/03/2022 02:31 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 546002

All characteristics are given for information only and may be changed

## STATION DE DIAGNOSTIC SPOT-CHECK 2VS PNI AVEC SOCLE A ROULETTES



Réf. 546002

### COMPOSITION

Code	Descriptive	Quantity
546000	Station De Diagnostic Spot-Check Vs Pni	1
546024	wheel stand for spot-check 2vs	1