

CUFF FOR ES 60/ADULT (M)

Réf. 550611



CE
0044



CUFF FOR ES 60/ADULT (M)



Réf. 550611

Grey nylon cuffs with quick coupling.

FEATURES			
Medical device	Yes	Classe I	
Accuracy measurement	Classe I		
Marquage	CE0044		
Barcode	3700446008932		
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	Width packaging (cm):		
Length (cm) :	Length packaging (cm) :		
Height (cm) :	Height packaging (cm) :		
Weight (kg) :	Weight packaging (kg) :		0,130
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)		
Temperature max	(N/A)		
Humidity min	(N/A)		
Humidity max	(N/A)		



26/04/2022 12:33 - This article is available on our website : www.gsh-med.fr with the following reference 550611

All characteristics are given for information only and may be changed