

## CUFF FOR ES 60/CHILD (S)

Réf. 550612



## CUFF FOR ES 60/CHILD (S)



**Réf. 550612**

Grey nylon cuffs with quick coupling.

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage	CE0044		Barcode	3700446008925
Guarantee (month):	6		Customs Code	9018901000
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :	0,090	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



26/04/2022 12:33 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 550612

All characteristics are given for information only and may be changed