

4VS MONITOR NIBP

Réf. 551600



CE
0123



4VS MONITOR NIBP



Réf. 551600

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446032944
Guarantee (month):	18	Customs Code	9018199560
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	3,760
SPECIFIC STORAGE CONDITION			
Temperature min	-40	Temperature max	55
Humidity min	5	Humidity max	95



31/03/2022 01:52 - This article is available on our website : www.gsh-med.fr with the following reference 551600

All characteristics are given for information only and may be changed