

SPO2 PEDIATRIC SENSOR FOR 4VS MONITOR



Réf. 551631



SPO2 PEDIATRIC SENSOR FOR 4VS MONITOR



Réf. 551631

FEATURES			
Medical device	No		
Marquage		Barcode	3700446031183
Guarantee (month):	6	Customs Code	90181990
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,093
SPECIFIC STORAGE CONDITION			
Temperature min	-40	Temperature max	55
Humidity min	5	Humidity max	95



31/03/2022 01:53 - This article is available on our website : www.gsh-med.fr with the following reference 551631

All characteristics are given for information only and may be changed