

## THERMAL PRINTER FOR 4VS SPENGLER MONITOR + PAPER ROLL



Réf. 551642



# THERMAL PRINTER FOR 4VS SPENGLER MONITOR + PAPER ROLL



Réf. **551642**

FEATURES			
Medical device	No		
Marquage		Barcode	3700446033118
Guarantee (month):	0	Customs Code	85049099
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,148
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:53 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 551642

All characteristics are given for information only and may be changed