

INFRARED THERMOMETER FOR 4VS MONITOR



Réf. 551644

INFRARED THERMOMETER FOR 4VS MONITOR



Réf. **551644**

FEATURES			
Medical device	No		
Marquage		Barcode	3700446033132
Guarantee (month):	0	Customs Code	9025192090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,162
SPECIFIC STORAGE CONDITION			
Temperature min	-40	Temperature max	55
Humidity min	5	Humidity max	95



31/03/2022 01:54 - This article is available on our website : www.gsh-med.fr with the following reference 551644