

## LIGHT BULB FOR CLAR MIRROR 3.5V (PREVIOUS MODEL 570101) SPENGLER

Réf. 570105



# LIGHT BULB FOR CLAR MIRROR 3.5V (PREVIOUS MODEL 570101) SPENGLER



**Réf. 570105**

Spare part for Clar mirror

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage			Barcode	3700446026936
Guarantee (month):	60		Customs Code	901890800009J
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :	0,005	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



31/03/2022 01:54 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 570105

All characteristics are given for information only and may be changed