

REPLACEMENT HEAD FOR OTOSCOPE SMARTLED 5500



Réf. 570518



REPLACEMENT HEAD FOR OTOSCOPE SMARTLED 5500



Réf. **570518**

FEATURES			
Medical device	No		
Marquage		Barcode	3700446026349
Guarantee (month):	0	Customs Code	90189042
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,046
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:54 - This article is available on our website : www.gsh-med.fr with the following reference 570518

All characteristics are given for information only and may be changed