

FETAL AND VASCULAR DOPPLER WITHOUT PROBE

Réf. 743301



FETAL AND VASCULAR DOPPLER WITHOUT PROBE



Réf. 743301

For prenatal ultrasound and examination of blood vessels. Choose the probes suited to your requirements. Large colour LCD screen displaying heart rate, battery indicator, pulse wave intensity and probe frequency. Supplied with two LR6 AA batteries. Ergonomic angled probes for a comfortable use.

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446004316
Guarantee (month):	24	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,430
SPECIFIC STORAGE CONDITION			
Temperature min	-10	Temperature max	55
Humidity min	5	Humidity max	93



26/10/2022 12:49 - This article is available on our website : www.gsh-med.fr with the following reference 743301

All characteristics are given for information only and may be changed