

PROBE FOR VASCULAR DOPPLER / 5 MHZ

Réf. 743311



PROBE FOR VASCULAR DOPPLER / 5 MHZ



Réf. 743311

5 MHz probe - vascular probe for doppler

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446008178
Guarantee (month):	12	Customs Code	90181990
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,145
SPECIFIC STORAGE CONDITION			
Temperature min	-10	Temperature max	55
Humidity min	5	Humidity max	93



31/03/2022 01:55 - This article is available on our website : www.gsh-med.fr with the following reference 743311

All characteristics are given for information only and may be changed