

PROBE FOR VASCULAR DOPPLER / 8 MHZ

Réf. 743312



PROBE FOR VASCULAR DOPPLER / 8 MHZ



Réf. 743312

8 MHz probe - cardiac probe for doppler

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446008161
Guarantee (month):	12	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,145
SPECIFIC STORAGE CONDITION			
Temperature min	-10	Temperature max	55
Humidity min	5	Humidity max	93



31/03/2022 01:55 - This article is available on our website : www.gsh-med.fr with the following reference 743312

All characteristics are given for information only and may be changed