

## ADULT TIPS FOR PEAK FLOWMETERS Ø 30 MM - BOX OF 100 PCS



Réf. L35587



# ADULT TIPS FOR PEAK FLOWMETERS Ø 30 MM - BOX OF 100 PCS



Réf. L35587

FEATURES			
Medical device	Yes	Classe I	
Marquage		Barcode	3700446001681
Guarantee (month):	0	Customs Code	48229000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,830
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:55 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference L35587

All characteristics are given for information only and may be changed