

## DISPLAY OF 12 THERMOMETRES TEMPO 10 BLUE

Réf. SP1623P



## DISPLAY OF 12 THERMOMETRES TEMPO 10 BLUE



Réf. **SP1623P**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446027513
Guarantee (month):	24	Customs Code	90251920
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,650
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:31 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference SP1623P

All characteristics are given for information only and may be changed

## DISPLAY OF 12 THERMOMETRES TEMPO 10 BLUE



**Réf. SP1623P**

### COMPOSITION

Code	Descriptive	Quantity
SP1620/	Thermometre Tempo 10 Blue	12