

DISPLAY OF 12 THERMOMETER TEMPO 10 FLEX GREEN

Réf. SP1644P



DISPLAY OF 12 THERMOMETER TEMPO 10 FLEX GREEN



Réf. **SP1644P**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446027537
Guarantee (month):	24	Customs Code	90251920
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	11,0
Length (cm) :		Length packaging (cm) :	14,0
Height (cm) :		Height packaging (cm) :	15,0
Weight (kg) :		Weight packaging (kg) :	0,624
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:32 - This article is available on our website : www.gsh-med.fr with the following reference SP1644P

All characteristics are given for information only and may be changed

DISPLAY OF 12 THERMOMETER TEMPO 10 FLEX GREEN



Réf. SP1644P

COMPOSITION

Code	Descriptive	Quantity
SP1643/	Thermometre Tempo 10 Flex Vert	12