

Réf. DSU230717



Réf. DSU230717

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage	CE0413		Barcode	3700446012304
Guarantee (month):	3	Customs Code		
Sterile	No			
INFORMATIONS	Main unit		Packaged unit	
Width (cm) :		Width packaging (cm):		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :		
SPECIFIC STORAGE CONDITION				
Temperature min		(N/A)	Temperature max	(N/A)
Humidity min		(N/A)	Humidity max	(N/A)



31/03/2022 01:56 - This article is available on our website : www.gsh-med.fr with the following reference DSU230717

All characteristics are given for information only and may be changed

2/2