

Réf. DSU230725



Réf. DSU230725

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE0413	Barcode	0110840935101000
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	15,5
Length (cm) :		Length packaging (cm) :	14,0
Height (cm) :		Height packaging (cm) :	3,5
Weight (kg) :		Weight packaging (kg) :	0,086
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:56 - This article is available on our website : www.gsh-med.fr with the following reference DSU230725