

OSCAR 2 POUCH

Réf. DSU230920

SunTech Medical®



CE



OSCAR 2 POUCH

Réf. **DSU230920**

Oscar 2TM other accessories. To carry the ABPM and for the patients' comfort.

FEATURES			
Medical device	Yes	Classe I	Accuracy measurement
Marquage	CE		Classe I
Guarantee (month):	0	Barcode	3700446027308
Sterile	No	Customs Code	90330000
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	2,0
Length (cm) :		Length packaging (cm) :	10,0
Height (cm) :		Height packaging (cm) :	12,0
Weight (kg) :		Weight packaging (kg) :	0,063
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:57 - This article is available on our website : www.gsh-med.fr with the following reference DSU230920

All characteristics are given for information only and may be changed