

Réf. DSU231000



## Réf. DSU231000

Reusable nylon cuff sold by 5. Single tube, screw-on coupling.

| FEATURES                   |        |                         |                      |
|----------------------------|--------|-------------------------|----------------------|
| Medical device             | Yes    | Classe I                | Accuracy measurement |
| Marquage                   | CE0413 |                         | Classe I             |
| Guarantee (month):         | 24     | Barcode                 | 0120840935100529     |
| Sterile                    | No     | Customs Code            |                      |
| INFORMATIONS               |        | Main unit               |                      |
|                            |        | Packaged unit           |                      |
| Width (cm) :               |        | Width packaging (cm):   |                      |
| Length (cm) :              |        | Length packaging (cm) : |                      |
| Height (cm) :              |        | Height packaging (cm) : |                      |
| Weight (kg) :              |        | Weight packaging (kg) : | 0,203                |
| SPECIFIC STORAGE CONDITION |        |                         |                      |
| Temperature min            | (N/A)  | Temperature max         | (N/A)                |
| Humidity min               | (N/A)  | Humidity max            | (N/A)                |



31/03/2022 01:57 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DSU231000

All characteristics are given for information only and may be changed