

Réf. DSU231001



Réf. DSU231001

Reusable nylon cuff sold by 5. Single tube, screw-on coupling.

FEATURES					
Medical device	Yes	Classe I	Accuracy measurement	Classe I	
Marquage	CE0413		Barcode	0120840935100383	
Guarantee (month):	24		Customs Code		
Sterile	No				
INFORMATIONS		Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):			
Length (cm) :		Length packaging (cm) :			
Height (cm) :		Height packaging (cm) :			
Weight (kg) :		Weight packaging (kg) :		0,216	
SPECIFIC STORAGE CONDITION					
Temperature min		(N/A)	Temperature max		(N/A)
Humidity min		(N/A)	Humidity max		(N/A)



31/03/2022 01:57 - This article is available on our website : www.gsh-med.fr with the following reference DSU231001

All characteristics are given for information only and may be changed