

Réf. DSU231004



## Réf. DSU231004

Reusable nylon cuff sold by 5. Single tube, screw-on coupling.

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage	CE0413		Barcode	0120840935104015
Guarantee (month):	24		Customs Code	
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	14,0	
Length (cm) :		Length packaging (cm) :	18,0	
Height (cm) :		Height packaging (cm) :	9,0	
Weight (kg) :		Weight packaging (kg) :	0,352	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)		Temperature max	(N/A)
Humidity min	(N/A)		Humidity max	(N/A)



31/03/2022 01:57 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DSU231004

All characteristics are given for information only and may be changed