

Réf. DSU231005



Réf. DSU231005

Reusable nylon cuff sold by 5. Single tube, screw-on coupling.

FEATURES			
Medical device	Yes	Classe I	Accuracy measurement
Marquage	CE0413		Classe I
Guarantee (month):	24	Barcode	0120840935104022
Sterile	No	Customs Code	
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	9,0
Length (cm) :		Length packaging (cm) :	18,0
Height (cm) :		Height packaging (cm) :	14,0
Weight (kg) :		Weight packaging (kg) :	0,450
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:57 - This article is available on our website : www.gsh-med.fr with the following reference DSU231005

All characteristics are given for information only and may be changed