

Réf. DSU231008



Réf. DSU231008

Reusable nylon cuff sold by 5. Single tube, screw-on coupling.

FEATURES			
Medical device	Yes	Classe I	Accuracy measurement
Marquage	CE0413		Classe I
Guarantee (month):	24	Barcode	3700446013226
Sterile	No	Customs Code	
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,770
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:58 - This article is available on our website : www.gsh-med.fr with the following reference DSU231008

All characteristics are given for information only and may be changed