

Réf. DSU231014



## Réf. DSU231014

Reusable nylon cuff sold by 5. Single tube, bayonet coupling.

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage	CE0413		Barcode	10840935104152
Guarantee (month):	24		Customs Code	
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :	0,333	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)		Temperature max	(N/A)
Humidity min	(N/A)		Humidity max	(N/A)



31/03/2022 01:58 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DSU231014

All characteristics are given for information only and may be changed