

Réf. DSU231025



Réf. DSU231025

Reusable nylon cuff sold by 5. Single tube, no coupling.

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage	CE0413		Barcode	3700446013738
Guarantee (month):	24		Customs Code	
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):	14,0	
Length (cm) :		Length packaging (cm) :	18,0	
Height (cm) :		Height packaging (cm) :	9,0	
Weight (kg) :		Weight packaging (kg) :	0,432	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



31/03/2022 01:58 - This article is available on our website : www.gsh-med.fr with the following reference DSU231025