

Réf. DSU231026



Réf. DSU231026

Reusable nylon cuff sold by 5. Single tube, no coupling.

FEATURES			
Medical device	No		
Marquage	CE0413	Barcode	0120840935104336
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,524
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:58 - This article is available on our website : www.gsh-med.fr with the following reference DSU231026

All characteristics are given for information only and may be changed