

Réf. DSU231030



Réf. DSU231030

| FEATURES | | | | |
|----------------------------|-----------|-------------------------|----------------------|---------------|
| Medical device | Yes | Classe I | Accuracy measurement | Classe I |
| Marquage | CE0413 | | Barcode | 3700446013820 |
| Guarantee (month): | 24 | | Customs Code | |
| Sterile | No | | | |
| INFORMATIONS | Main unit | Packaged unit | | |
| Width (cm) : | | Width packaging (cm): | | |
| Length (cm) : | | Length packaging (cm) : | | |
| Height (cm) : | | Height packaging (cm) : | | |
| Weight (kg) : | | Weight packaging (kg) : | | |
| SPECIFIC STORAGE CONDITION | | | | |
| Temperature min | (N/A) | Temperature max | (N/A) | |
| Humidity min | (N/A) | Humidity max | (N/A) | |



31/03/2022 01:59 - This article is available on our website : www.gsh-med.fr with the following reference DSU231030

All characteristics are given for information only and may be changed