

Réf. DSU231037



**Réf. DSU231037**

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage	CE0413		Barcode	0120840935104114
Guarantee (month):	24		Customs Code	
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :		
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



31/03/2022 01:59 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DSU231037

All characteristics are given for information only and may be changed