

## TIGHT 38-50CM X20

 SunTech Medical®

Réf. DSU231089



**CE**  
0413

DISPOSITIF MEDICAL  
MEDICAL DEVICE **CLASSE I**

## TIGHT 38-50CM X20

Réf. DSU231089

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE0413	Barcode	3700446014551
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	2,880
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:01 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DSU231089