

SPO2 REUSABLE DIGITAL SENSOR FOR ADULT 3 METERS

Réf. DSU250899

SunTech Medical®



SPO2 REUSABLE DIGITAL SENSOR FOR ADULT 3 METERS

Réf. **DSU250899**

FEATURES					
Medical device	Yes	Classe I	Accuracy measurement	Classe I	
Marquage	CE0413		Barcode	0108698870322340	
Guarantee (month):	3		Customs Code	90181990	
Sterile	No				
INFORMATIONS		Main unit	Packaged unit		
Width (cm) :			Width packaging (cm):	16,0	
Length (cm) :			Length packaging (cm) :	13,0	
Height (cm) :			Height packaging (cm) :	3,0	
Weight (kg) :			Weight packaging (kg) :	0,115	
SPECIFIC STORAGE CONDITION					
Temperature min		(N/A)	Temperature max		(N/A)
Humidity min		(N/A)	Humidity max		(N/A)



31/03/2022 02:02 - This article is available on our website : www.gsh-med.fr with the following reference DSU250899