

Réf. DSU551401



**Réf. DSU551401**

<BODY bgColor=#ffffff></BODY>

FEATURES					
Medical device	Yes	Classe IIa	Accuracy measurement	Classe II	
Marquage	CE0413		Barcode	0110840935103776	
Guarantee (month):	36		Customs Code	9018199560	
Sterile	No				
INFORMATIONS		Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):		23,5	
Length (cm) :		Length packaging (cm) :		19,5	
Height (cm) :		Height packaging (cm) :		14,8	
Weight (kg) :		Weight packaging (kg) :		0,500	
SPECIFIC STORAGE CONDITION					
Temperature min		(N/A)	Temperature max		(N/A)
Humidity min		(N/A)	Humidity max		(N/A)



31/03/2022 02:02 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DSU551401

All characteristics are given for information only and may be changed