

## MODULE TEMPERATURE COVIDIEN FILAC 3000 FOR STATION SUNTECH CT40

 SunTech Medical®

Réf. DSU551403



**CE**  
0413

DISPOSITIF MEDICAL CLASSE II  
MEDICAL DEVICE II

# MODULE TEMPERATURE COVIDIEN FILAC 3000 FOR STATION SUNTECH CT40

Réf. **DSU551403**

<BODY bgColor=#ffffff></BODY>

FEATURES				
Medical device	Yes	Classe IIa	Accuracy measurement	Classe II
Marquage	CE0413		Barcode	3700446032432
Guarantee (month):	36		Customs Code	9018.19.5500
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :			Width packaging (cm):	23,5
Length (cm) :			Length packaging (cm) :	19,5
Height (cm) :			Height packaging (cm) :	14,5
Weight (kg) :			Weight packaging (kg) :	0,612
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



02/12/2022 12:30 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DSU551403