

Réf. DSU551400



**Réf. DSU551400**

Using the SunTech technology for measuring blood pressure, the Pro VS Check® is the only scalable diagnostic station of this type. Designed to be easily used and adapted in all health care units, the modular design of this device allows you to add blood saturation and thermometry modules tailored to the needs of your patients. **RELIABILITY** The manual mode of measuring blood pressure allows reading and checking measurement in real-time, without removing the cuff. **PRECISION** The customizable average measurement mode of blood pressure ensures accurate and consistent readings. **CONNECTIVITY** Data transfer using IHE protocols with network connection via an integrated Ethernet port or optional WIFI. **FLEXIBILITY** The modular design allows you to add optional temperature and SpO2 modules. **SIMPLICITY** Intuitive user interface with LCD display.

FEATURES			
Medical device	Yes	Classe IIa	Accuracy measurement
Marquage	CE0413		Classe II
Guarantee (month):	36	Barcode	3700446028800
Sterile	No	Customs Code	9018905040
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	26,0
Length (cm) :		Length packaging (cm) :	31,0
Height (cm) :		Height packaging (cm) :	27,0
Weight (kg) :		Weight packaging (kg) :	3,200
SPECIFIC STORAGE CONDITION			
Temperature min		Temperature max	(N/A)
Humidity min		Humidity max	(N/A)



31/03/2022 01:29 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DSU551400

All characteristics are given for information only and may be changed