

Réf. DSU551455



Réf. DSU551455

FEATURES					
Medical device	Yes	Classe I	Accuracy measurement	Classe I	
Marquage	CE0413		Barcode	790069391859	
Guarantee (month):	3		Customs Code		
Sterile	No				
INFORMATIONS	Main unit		Packaged unit		
Width (cm) :		Width packaging (cm):			
Length (cm) :		Length packaging (cm) :			
Height (cm) :		Height packaging (cm) :			
Weight (kg) :		Weight packaging (kg) :		0,115	
SPECIFIC STORAGE CONDITION					
Temperature min		(N/A)		Temperature max	(N/A)
Humidity min		(N/A)		Humidity max	(N/A)



31/03/2022 02:02 - This article is available on our website : www.gsh-med.fr with the following reference DSU551455

All characteristics are given for information only and may be changed

2/2