

PLATE FOR PRINTER FOR TOURNIQUET HPM SERIE 800



Réf. GHP317009

PLATE FOR PRINTER FOR TOURNIQUET HPM SERIE 800



Réf. **GHP317009**

FEATURES			
Medical device	No		
Marquage		Barcode	3700446041304
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	19,0
Length (cm) :		Length packaging (cm) :	23,0
Height (cm) :		Height packaging (cm) :	18,5
Weight (kg) :		Weight packaging (kg) :	1,372
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:40 - This article is available on our website : www.gsh-med.fr with the following reference GHP317009

All characteristics are given for information only and may be changed