

## AUTOCLAVABLE CUFF ADULT LONG LEG 93X13CM (BLADDER 59X9 CM)

Réf. GHP552807

**Hpm**<sup>TM</sup>



# AUTOCLAVABLE CUFF ADULT LONG LEG 93X13CM (BLADDER 59X9 CM)



Réf. **GHP552807**

Nylon cuff with single silicone bladder suitable for autoclaving. With male Luer lock connector and cap for cleaning and sterilisation.

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | No        |                         |               |
| Marquage                   | CE        | Barcode                 | 3700446033712 |
| Guarantee (month):         | 12        | Customs Code            |               |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               |           | Width packaging (cm):   |               |
| Length (cm) :              |           | Length packaging (cm) : |               |
| Height (cm) :              |           | Height packaging (cm) : |               |
| Weight (kg) :              |           | Weight packaging (kg) : | 0,500         |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | (N/A)     | Temperature max         | (N/A)         |
| Humidity min               | (N/A)     | Humidity max            | (N/A)         |



31/03/2022 02:05 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference GHP552807

All characteristics are given for information only and may be changed