

AUTOCLAVABLE CUFF ADULT EXTRA LONG LEG 102X13CM(BLADDER 84X9CM)

Réf. GHP552808

Hpm™



AUTOCLAVABLE CUFF ADULT EXTRA LONG LEG 102X13CM(BLADDER 84X9CM)



Réf. **GHP552808**

Nylon cuff with single silicone bladder suitable for autoclaving. With male Luer lock connector and cap for cleaning and sterilisation.

FEATURES			
Medical device	No		
Marquage	CE	Barcode	3700446033729
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,630
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:05 - This article is available on our website : www.gsh-med.fr with the following reference GHP552808